



Protein Purification with the Äkta Avant

© Michael Jeltsch & University of Helsinki

This guide describes how to use the Äkta Avant 25, which is available for protein purification and located in the Protein Production and Purification Core Facility (<http://research.med.helsinki.fi/corefacilities/akta/index.html>). This manual is a work in progress and if you find anything missing, please let us know. This manual does not go into theory at all, but only describes the steps needed to perform a two step purification. In this guide, two often used examples of affinity chromatography methods are described: Ni²⁺NTA affinity chromatography and Protein A affinity chromatography. In our examples, we purify secreted proteins from conditioned cell culture supernatant. Because we have a rather large sample volume (400 ml), the sample is loaded directly to the column with the sample pump in the first (affinity chromatography) step. The second step is a size exclusion chromatography/gel filtration. In this step, the sample (the peak from the elution of the first chromatography step) is loaded with a syringe into the capillary loop, from which it is injected into the gel filtration column.

1. Buffers needed

For Ni²⁺NTA affinity chromatography

- Ni²⁺NTA binding buffer: 30 mM phosphate buffer pH 8.0, 400 mM NaCl
- Ni²⁺NTA elution buffer: 30 mM phosphate buffer pH 8.0, 400 mM NaCl, 250 mM imidazole
- column storage buffer: 20% ethanol

For Protein A affinity chromatography

- PBS-T (0.1% Tween-20)
- PAS elution buffer: 0.1 M citrate pH 3.0
- column storage buffer: 20% ethanol

2. Columns needed

For Ni²⁺NTA affinity chromatography

Depending on the amount of protein that you want to purify:

HisTrap FF 1 ml (GE Healthcare, 17-5319-01) for up to 40 mg of protein; HisTrap FF 5 ml (GE Healthcare, 17-5255-01) for more than 40 mg. Alternative columns for special cases are: HisTrap FF Crude (for purification from cell lysate), HisTrap HP (higher binding capacity at same column dimensions), HisTrap Excel (chemically more stable and higher affinity for his-tagged proteins, but only about 25% binding capacity)

For Protein A affinity chromatography

Depending on the amount of protein that you want to purify:

HiTrap rProtein A FF 1 ml (GE Healthcare, 17-5079-01) for up to 40 mg of human IgG₁;
HiTrap rProtein A FF 5 ml (GE Healthcare, 17-5080-01) for more than 40 mg of human IgG₁

3. Reservation, logging in and starting the program

3.1 Do not forget to reserve the machine at <http://jeltsch.akta.org/> → 'Online reservation



system'. You need to log in with your University login/password in order to make a reservation. Please leave your phone number in case we need to contact you.

- 3.2 Switch on the Äkta Avant and the computer if they are not running already. It does not matter in which order you switch them on. The Äkta ON/OFF switch is on the left side of the machine, near the bottom edge:
- 3.3 You need to login twice: First into the computer and then into the UNICORN program. To login into the computer, you may use either your own university account or if more than one person need to control the FPLC device, you can login in using a local account (user: Unicorn, pw: see the note attached to the computer). Then start the UNICORN program via the desktop icon. To login to the program, you need to use your own university account. In the login window, make sure that the Administrator box is not checked:
- 3.4 The machine should now already connect to the computer. If this doesn't happen automatically (usually the first time you log in) you can connect manually by clicking 'system' → 'Connect to systems' from the Unicorn program.

4. Preparation of the system for a run

- 4.1 Rinse the inlet tubes by dipping them in clean PBS/H₂O before placing them in your buffer bottles. Put your buffers in their appropriate places and inlet tubes in the right buffer bottles.

	Ni²⁺NTA affinity chromatography	Protein A affinity chromatography
A1:	Ni ²⁺ NTA binding buffer	PBS-T
B1:	Ni ²⁺ NTA elution buffer	0.1 M citrate pH 3.0
B3:	column storage buffer	column storage buffer
Buff:	Ni ²⁺ NTA binding buffer	PBS-T
S4:	sample	sample


- 4.2 Place the sample bottle into an ice-filled styrofoam box on top of the Äkta. You can additionally add 0.1% NaN₃ to the sample in order to prevent microbial growth.
- 4.3 Place a bottle into the fridge to collect your flowthrough, connect tubing Out5 to the bottle.

5. Priming/purging the pumps

- 5.1 If the inlet tubings A1, B1, B3 and Buff are NOT filled with buffer, you MUST fill them! Otherwise the pump will run dry and break! To do this you attach a syringe to all 6 pump purge valves and draw about 20 ml buffer.
In order to fill inlet tubing B3 and S4, you need to switch from B1 to B3 and from Buf to S3, respectively, using the UNICORN software. For B3: Manual instructions → Flow path → Pump B inlet: choose inlet B3. For S4: Manual instructions → Flow path → Sample inlet: choose inlet S4.
- 5.2 Inlet tubings A1, B1 and Buff are the default positions of the valves and can be filled without changing valve positions with the UNICORN software. The sample you draw into your syringe while filling tube S4 can be returned into the sample bottle.

6. Attaching the column

- 6.1 Start a slow flow (0.2 - 0.5ml/min) with the same buffer that has been used for column



storage (mostly 20% EtOH). Select the column position (1 to 5) that you want to use for connecting the column (the ready default methods in this guide all use position 1). Before attaching the new column fill the top of the column by slowly dripping buffer on it to remove any air that might be trapped inside.

6.2 Insert and screw in the top ferrule avoiding to insert any air.

6.3 Insert and screw in the top ferrule into the bottom of the column.

7. Calibrate the pH monitor

Calibrating the pH monitor is not necessary if you just want to monitor the pH, but don't need absolute pH values. We do it about once a month, but this means that it may be off by 0.1 to 0.4 pH units if you leave it uncalibrated.

8. Preparing the fraction collector

Fractions can be collected into many different tube types. The fraction collector holds six cassettes and we have cassette types available for the following tubes/plates:

- 4 cassettes for 96-deep-well plates
- 2 cassettes for 5-ml-tubes (40 tubes/cassette)
- 2 cassettes for 15-ml-tubes (15 tubes/cassette)
- 2 cassettes for 50-ml-tubes (6 tubes/cassette)


For the Ni²⁺+NTA affinity chromatography, we are using 15 ml Falcon tubes and for the Protein A affinity chromatography, we use Nunc 5 ml 75x12 Minisorp tubes. For the gel filtration, we use 96-deep-well plates with 1-ml-wells or 2-ml-wells. For a typical affinity purification you should place at least 10 tubes into the cassettes starting from position A1. For the Protein A affinity chromatography, you need to pipette 200 µl 1M Tris pH 8.5 into each fraction tube, because you need to neutralize the low pH of the elution buffer as soon as possible after the elution in order to prevent damaging your protein. For the gel filtration, place one 96-deep-well block into a 96-plate cassette.

9. Creating a method

A description of method creation is coming in a later version of this manual. For routine purifications (including both Ni²⁺+NTA and Protein A affinity chromatography), there are ready methods available in the Unicorn software, which can be used: When you start the run, choose methods Ni²⁺+NTA_HisTrap_1ml_default or Ni²⁺+NTA_HisTrap_5ml_default for the Ni²⁺+NTA affinity chromatography depending on the column you are using and HiTrapPAS_1ml_default or HiTrapPAS_5ml_default for the Protein A affinity chromatography depending on the column you are using. For a standard gel filtration on the Superdex 200 Increase column with a 0.5-ml-capillary loop, you should use the method *SEC200Increase_default_no_eqilibration* if the column has been already equilibrated in running buffer and *SEC200Increase_default_with_eqilibration* if the column is still in storage buffer (20% ethanol).

10. Sample loading into the capillary loop (only for step 2: gel filtration)

For the gel filtration run, you need to manually load the capillary loop with your sample before you start the run. Depending on your sample volume, you need to choose a volume-matched capillary loop. If you want to maximize your recovery, always use a capillary loop that is at least 30% larger than the volume of your sample. However, if peak separation is your goal, use a very small sample volume, a very small capillary loop and overfill it (e.g. fill a 0.5-ml-loop with 0.8 ml sample and resolve it on the HiLoad 16/600 Superdex pg).

- 
- 10.1 Attach the correct capillary loop to positions LoopF and LoopE. Wash the loop by connecting a syringe to position Syr and pushing about 10 loop volumes running buffer through the loop.
 - 10.2 Load your sample into a syringe (e.g. a 1-ml-syringe) avoiding/removing any air bubbles and insert the syringe to position Syr
 - 10.3 Push the syringe content SLOWLY into the sample loop. Excess sample exits at position W1 and you may recover it.

11. Starting the method run

- 11.1 Click the “folder” in the top left corner of the program window and choose a ready method.
- 11.2 When the program asks for the column information use the bar code scanner to scan the bar code on the column or find it on the list.
- 11.3 Double check that all tubes are in correct buffers.
- 11.4 Start the run by clicking the 'play' button.
- 11.5 After the run has finished, you can remove your column. The default methods equilibrate the column in the end to 20% ethanol for storage. Even though it is sold as a single-use device, you can reuse it a couple of times. However, if you use it for a different protein, you need to do rigorous cleaning in between. Deattach the column from its tubings and firmly screw stop plugs to both ends to prevent the column from drying out. Store the column at +4°C.

12. Collecting, pooling and concentrating fractions from the affinity chromatography (first step)

After the run has finished, you can look up from the evaluation module, in which fraction(s) the protein eluted and how concentrated it is. You can pool the desired fractions. If your pool has a large volume, you need to concentrate it before you load it to the gel filtration column. Typical gel exclusion columns (Superdex 10/300 GL or Superdex Increase) accept a maximum sample volume of ~0.5 ml (for volumes of ~0.5 to ~1.5 ml, we have the HiLoad 16/600 Superdex pg and for volumes of ~1.5 to ~5 ml, we have the HiLoad 26/600 Superdex pg). Depending on the size of your protein, you should use either the Superdex 75 (proteins from 3 to 70 kDa) or the Superdex 200 (proteins from 10 to 600 kDa). Concentration is done most conveniently by ultrafiltration (e.g Amicon Ultra-4 10K for up to 4 ml, Amicon Ultra-15 10K for up to 15 ml). Because of the high salt and imidazole concentration, it is recommended to dilute the pooled elution fractions 1:1 with water before concentration. Be careful not to overconcentrate (because your protein might precipitate). Overconcentration can happen locally near to the ultrafiltration membrane if you try to concentrate in one centrifugation step. It is better to do several shorter spins and homogenize the concentrate in between to avoid high protein concentrations near the membrane.