

Sekvensointikortti ABI-310

dsDNA 200-500 ng
ssDNA 50-100 ng

PCR products:

100-200 bp 5 ng
200-500 bp 10 ng
500-1000 bp 20 ng
2000-3000 bp 50 ng

Pvm: _____

Tilaaja: _____

Primer 3.2 pmol Total vol 12 ul in aqua.

Samples in 0.5 ul PCR tubes.

	Template	ul	Primer
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			